



Before using, position the eFOLDi Lite in an open area away from pets and small children. Unfolding is easy, and with practice, you can set it up in seconds. We recommend practicing before your first outing.



HOW TO UNFOLD

STEP 1 - RELEASE THE STEERING COLUMN

Remove all packaging before you begin.



Push the yellow latch to release the steering column, allowing it to swing free.

Lift the handlebars outward and upward until they click securely into a vertical position.

STEP 2 - UNLOCK & UNFOLD



Hold the steering column & locate the small yellow toggle on the side of the Lite and pull it to unlock the front wheel assembly.

Hold the steering column and seat back, and pull them apart. The scooter will unfold and you will hear a click as it locks into its open position.

STEP 3 - FINAL SAFETY CHECKS



Ensure the locking latches are securely fastened to the steel bar by pressing your foot down on the footplate.



Check that the red magnetic brake lever is in the LOCKED position to engage the motor.

GETTING STARTED



To turn on the scooter, press the grey button below the battery screen. Select (D) for Drive or (R) for Reverse, then press the thumb throttle to go.

Be sure to re-charge the battery when the display shows amber.



To turn the scooter into rolling luggage pull the handle bar up and press the yellow lever, securing it on the second stop point..



To charge, connect the power supply cable to the charging port located at the top of the steering column.



To remove the battery lift the latch and carefully pull it upwards.

The charging port is located on the side of the battery.

HOW TO FOLD TO STORE

STEP 1 - FOLD THE BACK REST & STEERING COLUMN



Fold the seat backrest down flat against the seat base.



Press and hold the small yellow latch on the front of the steering column, then fold it down until the handle bar down until they click into place, secured by the magnet.



With armrests fitted, pull the black toggle to release the seat and fold.

STEP 2 - COLLAPSE & LOCK THE SCOOTER



Grip the scooter by the steering column and carry handle & lift the seat frame to an upright position.



Tilt the scooter back at an angle (you can use your foot behind the wheel to steady it if needed), balancing it onto the casters.



Press the steering column inward until the front wheel and latch click securely into place. Then, push the column firmly to lock it tightly against the frame.



Please note this leaflet is intended to be a quick guide only. For full instructions, please refer to our user product manual. If you have any questions, please call us on 0800 915 0800 or email services@efoldi.com