

Explorer



Rediscover your freedom...

Ultra-lightweight, folding lifestyle scooter



Scan for more information



eFOLDi



...and travel with confidence.

The eFOLDi Explorer is the ultimate solution for inclusive independence

As a recognised leader in the design and manufacture of portable scooters and powered wheelchairs, eFOLDi strives to help people rediscover their freedom and independence.

eFOLDi's award-winning mobility solutions are not only easy to store and transport, they help make travel simple and achievable – whether that's by air, land or sea.



eFOLDi



Fold and go further by car, bus, train, plane, or cruise liner – it's your choice

The award-winning eFOLDi Explorer is an ultra-portable lifestyle scooter for road* or pavement use.

Weighing just 17kg**, it folds and unfolds in a matter of seconds, so is easy to store and transport.

No complex assembly is required – simply, unfold and go.

What makes eFOLDi unique...

Lightweight – carry and lift with ease



Compact – folds down to cabin luggage size



Portable – easy to take with you, anywhere



Robust – durable materials give peace of mind



Explorer

Compact, comfortable, quick and simple to fold...the Explorer combines it all

- Easy to use, folding lifestyle scooter, for road* or pavement use
- Exceptionally lightweight at only 17kg**, yet remarkably robust carrying passengers up to 120kg/19 Stone
- Nimble and very manoeuvrable
- Folds with ease, fits into most car boots and is Airline-safe
- Travels up to 8mph with a range of 14 miles
- The perfect travel companion!



The lighter way to empower your life

The perfect blend of uncompromising style, portability and versatility

"It's simply magic. My eFOLDi Explorer gives me the complete independence to travel and do all the things I used to be able to do when I was younger and more mobile."

Maggie, London





Robust, carries users up to 120kg/265lbs

Ignition key for added security and peace of mind

Supportive cushion with suspension

Simple, responsive controls

Pneumatic tyres for a smooth ride

Powerful motor for varying terrain



Spacious storage bag

Easy-grip pull handle

Easy to lift, store and pull-along

Airline-friendly, lithium battery

Highly visible brake and indicator lights

Folds into a compact size within seconds



Specifications*

Maximum speed	8mph / 12kph
Maximum range	14 miles / 22km
Maximum user weight	19 Stone / 120kg / 265lbs
Climbing ability	1:6 gradient
Folded dimensions	26 x 22 x 16 inches / 66 x 56 x 40cm
Weight (excl. battery)	37.5lbs / 17kg
Weight (incl. battery)	42lbs / 19kg
Type of use	Road (UK legal) / pavement
Braking system	Manual
Airline friendly	Yes

* Specifications may vary due to manufacturing tolerances

Explorer

...and enjoy being you.

"My eFOLDi Explorer is an absolute lifesaver – it makes life so easy for me."

Danielle, Colchester



Scan for more information

0800 999 8123

sales@efoldi.com www.efoldi.com

Your local eFOLDi retailer:



SunTech UK Limited t/a eFOLDi, 25 Ormside Way, Holmethorpe Industrial Estate, Redhill, Surrey RH1 2LW.

eFOLDi