



unfold the future!



**Your best travel companion
wherever you go!**

Foldable frame

**Telescopic and
turnover footrest**

**Detachable
lithium battery**

**Quick release for
backrest folding**

**Controller with
rotary joint**



SPECIFICATION

UNFOLDED DIMENSIONS(MM)

Motor	2 x 250W brush motor	LENGTH	WIDTH	HEIGHT
Controller	Micon			
Drive Type	RWD			
Battery	24V 10AH x 2PC Lithium battery	1070	680	950
Charger	DC24V/2.0A	FOLDED DIMENSIONS(MM)		
Feame	Aluminum with anodizing	LENGTH	WIDTH	HEIGHT
Front Wheel	Solid PU wheel in 2*8 and Nylon hub			
Rear Wheel	Solid PU in 2*12.5 and Aluminum hub with rubber tyre			
Brake Tybe	Magnetic brake	800	680	460
Travel Distance	≥15KM	NET/GROSS WEIGHT(KG)		
Max Speed	6KM/H	PACKING SIZE(MM)		
Weight Capacity	125KG	25	28	820X780X400
Spanning Ability	100MM			
Turning Radius	900MM			
Maximum Climbing Capacity	10°±1°			

www.efoldi.co.uk

eFOLDi

Folding Power Chair



FEATURES

- Lightweight and foldable for easy handling, storage and transportation
- High strength folding aluminum frame with anodizing
- Suspension system that absorbs impact and provides comfort and stability
- Rotatable and telescopic armrests and controller lock lever for comfortable use
- Telescopic and turnover footrest
- Adjustable backrest angle in three positions
- Optional pneumatic tyre and nylon hub

More award winning eFOLDi products



Tel: +44(0)20 3143 5168 (London)
+44(0)13 3343 9356 (Scotland)

Email: hello@efoldi.com
www.efoldi.co.uk