

Lightweight mobility solutions

Rediscover your freedom





Recognised leaders in the design and manufacture of ultra-portable scooters and power chairs

Explorer

My name is Sumi, and welcome to eFOLDi - Home to some of the world's lightest and most compact folding mobility scooters and power chair.

We strive to help people rediscover their freedom and independence.

Lite

When my father lost his mobility in a work-place accident some years ago, he simply couldn't find a mobility scooter anywhere, that was easy to use, safe and comfortable - but also lightweight and portable. Undeterred, my father set about designing his own scooter - and that is how the first eFOLDi mobility scooter came to be born!

Power Chair

Since then, eFOLDi has developed into a thriving family business, designing and manufacturing innovative mobility products that help make travel simple, achievable and most importantly, fun! The eFOLDi scooter was also a proud winner of the prestigious 'British Invention Of The Year' award in 2016. Legendary innovator and entrepreneur, Sir Richard Branson, even named the eFOLDi scooter "folding magic".

Join us in making travel easy and enjoyable again, and let an eFOLDi be your perfect travel companion.

As seen on TV

Sumi Wang

Founder



What makes eFOLDi unique



Lightweight

Carry and lift with ease



Portable

Easy to take with you, anywhere



Compact

Folds down to cabin luggage size



Robust

Peace of mind with strong, durable materials



“eFOLDi is folding magic.”

Sir Richard Branson





Explorer

BRITISH
INVENTION
OF THE YEAR
AWARD

Ultra-lightweight folding scooter

The multi-award-winning Explorer from eFOLDi is one of the world's lightest, most compact folding UK road-legal scooters. This ultra-portable scooter folds and unfolds simply, in a matter of seconds and requires no complex assembly - just unfold and go.

Lightweight • Portable • Compact • Robust



Lightweight
at only 17kg/37.5lbs*



Easy to use
controls



Compact size
folds with ease



Airline friendly battery
14 miles/22km
per charge



Travel in comfort
suspension/
pneumatic tyres



Fully UK road legal
8mph/12kmh
top speed

*excl. battery

“

My eFOLDi has impacted hugely on my life. It has enabled me to go out so much more than before. Now I have got my freedom and independence back, I can go out whenever and wherever I want.

Ray, London

”



Lite

NAIDEX
INNOVATION
OF THE YEAR
AWARD

Ultra-lightweight folding scooter

The eFOLDi Lite is one of the world's lightest, most compact folding travel-scooters. The Lite is the latest, and most advanced scooter from eFOLDi - weighing just 15kg*. It folds and unfolds in just a few seconds and is easy to store and transport. No complex assembly is required - just unfold and go.

Lightweight • Portable • Compact • Robust



Lightweight
at only 15kg/33lbs*



Easy to use
thumb controls



Compact size
folds with ease



Travel in comfort
suspension/
pneumatic tyres



Airline friendly battery
10 miles/16km
per charge



Robust
and carries up to
19 Stone/120kg/265lbs

*excl. battery



eFOLDi has changed my life so much for the better. It allows me to participate so much more in the world. It folds up extremely neatly like a small suitcase, and is light enough for someone to easily lift into a car boot. I have gone from being largely isolated to being much more able to participate in outings.

Jackie, Esher



Power Chair



Ultra-lightweight folding power chair

The eFOLDi Power Chair is one of the world's lightest, most compact folding powered wheelchairs. This ultra-portable, lightweight power chair is the perfect travel companion, and folds and unfolds in a matter of seconds. No complex assembly is required - just unfold and go.

Lightweight • Portable • Compact • Robust



Handy controls
for attendant to drive



Easy to use
joystick and control pad



Compact size
folds with ease



Travel in comfort
durable, breathable fabric



Airline friendly batteries
8 miles/13km per charge



Lightweight
at only 15kg/33lbs*

*excl. battery



With an eFOLDi you can just fold it up, and take it with you on buses, trains, or planes. Perfect for me to get out and about with my family and even go on holiday abroad!

Anna, Reigate



Specifications

Explorer



Lite



Power Chair



Maximum Speed	8mph / 12kph	4mph / 6kph	4mph / 6kph
Maximum Range	14 miles / 22km	10 miles / 16km	8 miles / 13km
Maximum User Weight	19 Stone / 120kg / 265lbs	19 Stone / 120kg / 265lbs	21 Stone / 135kg / 297lbs
Climbing Ability	1:6 Gradient	1:10 Gradient	1:10 Gradient
Folded Dimensions	26 x 22 x 16 inches 66 x 56 x 40cm	26 x 22 x 15 inches 65 x 55 x 37cm	22 x 13 x 36 inches 58 x 34 x 92cm
Weight (excl. battery)	37.5 lbs /17kg	33lbs /15kg	33lbs /15kg
Weight (incl. battery)	42lbs/19kg	37.5lbs/17kg	37.5lbs/17kg
Type of Use	Road (UK legal)/Pavement	Pavement Only	Pavement Only
Braking System	Manual	Electro-Magnetic	Electro-Magnetic
Airline friendly	Yes	Yes	Yes

Specifications may vary due to manufacturing tolerances.

Accessories

To complement your eFOLDi, we offer a wide range of accessories tailored to suit a variety of needs. This includes a host of storage solutions, protective covers, cushions and battery charger options.




Innovation

Recognised as leaders in the design and manufacture of ultra-portable scooters and power chairs, we strive to help people rediscover their freedom and independence.

Our products are manufactured using the highest quality materials to ensure they're robust enough to withstand the toughest of conditions and terrain. Our scooter frames are made from high-grade aluminium or magnesium alloys, our high-specification brushless motors deliver exceptional

performance, and our batteries use advanced lithium-ion technology, so they're not only durable, they offer outstanding levels of reliability.



Support

All eFOLDi products come with a comprehensive 12-month UK manufacturer's warranty. In the unlikely event of a problem, our UK support centre and servicing team are always on hand to help. They are trained to diagnose and rectify any problems, including warranty repairs and any ongoing service and maintenance requirements.



For more information...

Call us today on:
0800 999 8123

Or email us at:
sales@efoldi.com



Authorised eFOLDi Dealer:

eFOLDi 

0800 999 8123 • sales@efoldi.com • www.efoldi.com

SunTech UK Limited t/a eFOLDi, 25 Ormside Way, Holmethorpe Industrial Estate, Redhill, Surrey RH1 2LW.

UK
CA



FREEDOM • FRIENDSHIP • FUN