



7 Year
Battery + 2
Year Scooter
Warranty**

Whether it's the shops, the park, or a quick getaway — the Compass is ready to go.

Our first 4-wheel, two-step folding scooter, designed to make your journeys easier and more enjoyable.

Weighing just 11kg (without the battery), it's lightweight and perfect for confidently driving across different terrains.

Plus, it folds in seconds, making it simple to transport and store wherever you go.



Lightweight
at only
11kg/25lbs*




Top speed
4mph/6kph



Range
up to
10miles/16km



Warranty
Backed by
full support &
warranty


Ideal For: Those with minor mobility issues who need extra support and stability, seeking a scooter that offers increased confidence, security, and ease of use for getting around.



The lighter way to empower your life

The perfect blend of uncompromising style, portability and versatility



Flip-up armrests (optional)
for easy transfer in and out of the scooter

Ergonomic, breathable seating
with pressure-relieving memory foam

Automatic braking
slows/stops upon throttle release

Charging port
for on-or-off-board charging

Robust
carries users up to 120kg/265lbs/ 19 Stone

Bi-folding convenience
perfect for users with limited dexterity or strength

Powerful motor
low-vibration motor with whisper-quiet performance

Secure-Turn® technology
built-in gyroscope helps with stability and balance by detecting changes in orientation and movement



Simple & light-touch controls
keyless operation with easy one-touch button start

Flexible seat back
designed to accommodate different body shapes



Spacious foot area
Helps prevent cramping & discomfort

Automotive grade Lithium-ion battery
with 10 mile / 16km range, removable for on-or-off-board charging



Solid tyres
with no pumping requirement & multi-terrain suitability

Lightweight & portable
weighing just 11kg*



"It's so clever, the scooter adjusts automatically when I turn, which makes driving feel safe and smooth."

Ros, Verified Buyer

The Compass brings together comfort, compact design, and effortless folding — all in one smart scooter.

- Stable, comfortable and ready for every journey
- Exceptionally lightweight at only 11kg*, yet remarkably robust, carrying passengers up to 120kg/19 Stone
- Nimble and very manoeuvrable with Secure-Turn® technology for extra confidence and stability
- Folds in seconds to fit easily into car boots and is airline-friendly
- Travels up to 4 miles (6kph) an hour with a range of 10 miles (16km)
- The perfect travel companion!

Accessories

Make every journey easier with thoughtful extras designed for your comfort and convenience - from clever storage and protective covers to cushions, attachments, spare batteries and chargers.



*without battery/accessories
**terms and conditions apply

Wherever life takes you,
we'll be right here, making
sure you have the mobility
and freedom to enjoy
every moment.

As a recognised leader in the design and
manufacture of portable scooters and powered
wheelchairs, eFOLDi strives to help people
rediscover their freedom and independence.

eFOLDi's award-winning mobility solutions are not only easy
to store and transport, they help make travel simple and
achievable – whether that's by air, land or sea.



Compass Specifications*

Maximum Speed	4mph / 6kph
Maximum Range	10 miles / 16km
Maximum User Weight	19 Stone / 120kg / 265lbs
Climbing Ability	1:10 Gradient
Turning Radius	135cm (53*)
Folded Dimensions	92 x 43 x 40cm 36 x 17 x 16 inches
Unfolded Dimensions	92 x 43 x 96cm 36 x 17 x 38 inches
Weight (excl. battery)	11kg/25lbs
Weight (incl. battery)	13kg/29lbs
Type of Use	Pavement and indoor
Braking System	Electro-magnetic
Airline friendly	Yes

*Specifications may vary due to manufacturing tolerances

Why eFOLDi stands apart



Lightweight
Carry and lift
with ease



Warranty
Backed by full
support & warranty



Compact
Folds down to 'hold'
luggage size



Robust
Peace of mind with
strong, durable
materials



Innovation

Our products are manufactured using the highest
quality materials to ensure they're robust enough
to withstand the toughest of conditions and terrain.
Our frames are made from high-grade aluminium
or magnesium alloys, our high-specification
brushless motors deliver exceptional performance,
and our batteries use advanced lithium-ion
technology, so they're not only durable, they offer
outstanding levels of reliability.



7 Year Battery + 2 Year Scooter Warranty**

**terms and conditions apply



Scan for more
information



For more information...

Call us today on
0800 999 8123
Or email us at:
sales@efoldi.com



SunTech UK Limited t/a eFOLDi , 25 Ormside Way,
Holmehorpe Industrial Estate, Redhill, Surrey RH1 2LW.



eFOLDi Compass

Lightweight
mobility solutions

For wherever life takes you

www.efoldi.com



Scan for more
information