



eFOLDi Explorer

Lightweight mobility solutions

For wherever life takes you

www.efoldi.com



Scan for more
information




**7 Year
Battery + 2
Year Scooter
Warranty*****

Fold and go further by car, bus, train, plane, or cruise liner – it's your choice

The award-winning eFOLDi Explorer is an ultra-portable lifestyle scooter for road* or pavement use.

Weighing just 17kg**, it folds and unfolds in a matter of seconds, so is easy to store and transport.

No complex assembly is required – simply, unfold and go.



Ideal For: Individuals with minor mobility issues who want to travel further and faster, enjoying greater freedom and independence.



Lightweight
at only
17kg/38lbs**



Top speed
8mph/12kph



Range
up to
14miles/22km



Warranty
Backed by
full support &
warranty

*Approved for UK roads only
**without battery/accessories
***terms and conditions apply

The lighter way to empower your life

The perfect blend of uncompromising style, portability and versatility





"With the eFOLDi I'm able to visit places I haven't been able to get to for a while. I'm very happy to tell people – buy an eFOLDi."

Peter, Verified Buyer 

Compact, comfortable, quick and simple to fold... the Explorer combines it all.

- Easy to use, folding lifestyle scooter, for road* or pavement use
- Exceptionally lightweight at only 17kg**, yet remarkably robust carrying users up to 120kg/19 Stone
- Nimble and very manoeuvrable
- Folds with ease, fits into most car boots and is Airline-safe
- Travels up to 8mph (12kph) with a range of 14 miles (22km)
- The perfect travel companion!

Accessories

Make every journey easier with thoughtful extras designed for your comfort and convenience – from clever storage and protective covers to cushions, attachments, spare batteries and chargers.



**Wherever life takes you,
we'll be right here, making
sure you have the mobility
and freedom to enjoy
every moment.**

As a recognised leader in the design and manufacture of portable scooters and powered wheelchairs, eFOLDi strives to help people rediscover their freedom and independence.

eFOLDi's award-winning mobility solutions are not only easy to store and transport, they help make travel simple and achievable – whether that's by air, land or sea.

COMPASS



LITE



EXPLORER



NAVIGATOR



Why eFOLDi stands apart



Lightweight

Carry and lift
with ease




Warranty

Backed by full
support & warranty



Compact

Folds down to 'hold'
luggage size



Robust

Peace of mind with
strong, durable
materials



Innovation

Our products are manufactured using the highest quality materials to ensure they're robust enough to withstand the toughest of conditions and terrain. Our frames are made from high-grade aluminium or magnesium alloys, our high-specification brushless motors deliver exceptional performance, and our batteries use advanced lithium-ion technology, so they're not only durable, they offer outstanding levels of reliability.

Explorer Specifications*

Maximum Speed	8mph / 12kph
Maximum Range	14 miles / 22km
Maximum User Weight	19 Stone / 120kg / 265lbs
Climbing Ability	1:6 Gradient
Turning Radius	100cm (40*)
Folded Dimensions	66 x 56 x 40cm 26 x 22 x 16 inches
Unfolded Dimensions	121 x 56 x 83cm 48 x 22 x 33 inches
Weight (excl. battery)	17kg/38lbs
Weight (incl. battery)	19kg/42lbs
Type of Use	Road (UK legal)/Pavement, Indoor
Braking System	Mechanical Hand Brake
Airline friendly	Yes

*Specifications may vary due to manufacturing tolerances



7 Year Battery + 2 Year Scooter Warranty***

***terms and conditions apply



Scan for more
information

eFOLDi 

For more information...

Call us today on
0800 999 8123

Or email us at:
sales@efoldi.com

Your local eFOLDi retailer:



SunTech UK Limited t/a eFOLDi, 25 Ormside Way,
Holmethorpe Industrial Estate, Redhill, Surrey RH1 2LW.