



7 Year
Battery + 2
Year Scooter
Warranty**

Fold and go further by car, bus, train, plane, or cruise liner – it's your choice

The eFOLDi Lite is one of the world's lightest, most compact folding travel-scooters.

Weighing just 14kg*, the eFOLDi Lite folds and unfolds in just a few seconds so is easy to store and transport.

No complex assembly is required – simply, unfold and go.



Lightweight
at only
14kg/31lbs*



Top speed
4mph/6kph




Range
up to
10miles/16km



Warranty
Backed by
full support &
warranty

* without battery/accessories
**terms and conditions apply


Ideal For: Those with minor mobility issues who prioritise extra comfort and ease, looking for a scooter that supports their everyday travel needs with added convenience.



The lighter way to empower your life

The perfect blend of uncompromising style, portability and versatility



"What a fantastic device. It's been so amazing for my husband to have something so portable."

Lisa, Verified Buyer

Compact, comfortable, quick and simple to fold... the Lite combines it all.

- Easy to use, folding travel scooter
- Exceptionally lightweight at only 14kg*, yet remarkably robust carrying users up to 120kg/19 Stone
- Nimble and very manoeuvrable, with auto turning speed controller for extra stability
- Folds with ease, fits into most car boots and is Airline-safe
- The perfect travel companion!

Accessories

Make every journey easier with thoughtful extras designed for your comfort and convenience - from clever storage and protective covers to cushions, attachments, spare batteries and chargers.



Wherever life takes you,
we'll be right here, making
sure you have the mobility
and freedom to enjoy
every moment.

As a recognised leader in the design and
manufacture of portable scooters and powered
wheelchairs, eFOLDi strives to help people
rediscover their freedom and independence.

eFOLDi's award-winning mobility solutions are not only easy
to store and transport, they help make travel simple and
achievable – whether that's by air, land or sea.



Lite Specifications*

Maximum Speed	4mph / 6kph
Maximum Range	10 miles / 16km
Maximum User Weight	19 Stone / 120kg / 265lbs
Climbing Ability	1:10 Gradient
Turning Radius	100cm (40°)
Folded Dimensions	70 x 55 x 37cm 27.5 x 22 x 15 inches
Unfolded Dimensions	117 x 55 x 94cm 46 x 22 x 37 inches
Weight (excl. battery)	14kg/31lbs
Weight (incl. battery)	16kg/35lbs
Type of Use	Pavement and indoor
Braking System	Electro-magnetic
Airline friendly	Yes

*Specifications may vary due to manufacturing tolerances

Why eFOLDi stands apart



Lightweight
Carry and lift
with ease



Warranty
Backed by full
support & warranty



Compact
Folds down to 'hold'
luggage size



Robust
Peace of mind with
strong, durable
materials



Innovation

Our products are manufactured using the highest
quality materials to ensure they're robust enough
to withstand the toughest of conditions and terrain.
Our frames are made from high-grade aluminium
or magnesium alloys, our high-specification
brushless motors deliver exceptional performance,
and our batteries use advanced lithium-ion
technology, so they're not only durable, they offer
outstanding levels of reliability.



7 Year Battery + 2 Year Scooter Warranty**

**terms and conditions apply



Scan for more
information



For more information...

Call us today on
0800 999 8123
Or email us at:
sales@efoldi.com



SunTech UK Limited t/a eFOLDi, 25 Ormside Way,
Holmesthorpe Industrial Estate, Redhill, Surrey RH1 2LW.



eFOLDi Lite

Lightweight
mobility solutions

For wherever life takes you

www.efoldi.com



Scan for more
information