


# Lightweight mobility solutions

For wherever life takes you

[www.efoldi.com](http://www.efoldi.com)



Scan for more  
information





**7 Year  
Battery + 2  
Year Scooter  
Warranty\*\***

## Fold and go further by car, bus, train, plane, or cruise liner – it's your choice

The eFOLDi Lite is one of the world's lightest, most compact folding travel-scooters.

Weighing just 14kg\*, the eFOLDi Lite folds and unfolds in just a few seconds so is easy to store and transport.

No complex assembly is required – simply, unfold and go.



**Lightweight**  
at only  
14kg/31lbs\*



**Top speed**  
4mph/6kph



**Range**  
up to  
10miles/16km



**Warranty**  
Backed by  
full support &  
warranty

\* without battery/accessories

\*\*terms and conditions apply



**Ideal For:** Those with minor mobility issues who prioritise extra comfort and ease, looking for a scooter that supports their everyday travel needs with added convenience.



# The lighter way to empower your life

The perfect blend of uncompromising style, portability and versatility







"What a fantastic device.  
It's been so amazing for  
my husband to have  
something so portable."

Lisa, Verified Buyer 

## Compact, comfortable, quick and simple to fold... the Lite combines it all.

- Easy to use, folding travel scooter
- Exceptionally lightweight at only 14kg\*, yet remarkably robust carrying users up to 120kg/19 Stone
- Nimble and very manoeuvrable, with auto turning speed controller for extra stability
- Folds with ease, fits into most car boots and is Airline-safe
- The perfect travel companion!

## Accessories

Make every journey easier with thoughtful extras designed for your comfort and convenience - from clever storage and protective covers to cushions, attachments, spare batteries and chargers.





**Wherever life takes you,  
we'll be right here, making  
sure you have the mobility  
and freedom to enjoy  
every moment.**

As a recognised leader in the design and manufacture of portable scooters and powered wheelchairs, eFOLDi strives to help people rediscover their freedom and independence.

eFOLDi's award-winning mobility solutions are not only easy to store and transport, they help make travel simple and achievable – whether that's by air, land or sea.

COMPASS



LITE



EXPLORER



NAVIGATOR



## Why eFOLDi stands apart



### Lightweight

Carry and lift  
with ease




### Warranty

Backed by full  
support & warranty



### Compact

Folds down to 'hold'  
luggage size



### Robust

Peace of mind with  
strong, durable  
materials



## Innovation

Our products are manufactured using the highest quality materials to ensure they're robust enough to withstand the toughest of conditions and terrain. Our frames are made from high-grade aluminium or magnesium alloys, our high-specification brushless motors deliver exceptional performance, and our batteries use advanced lithium-ion technology, so they're not only durable, they offer outstanding levels of reliability.

## Lite Specifications\*

Maximum Speed	4mph / 6kph
Maximum Range	10 miles / 16km
Maximum User Weight	19 Stone / 120kg / 265lbs
Climbing Ability	1:10 Gradient
Turning Radius	100cm (40*)
Folded Dimensions	70 x 55 x 37cm 27.5 x 22 x 15 inches
Unfolded Dimensions	117 x 55 x 94cm 46 x 22 x 37 inches
Weight (excl. battery)	14kg/31lbs
Weight (incl. battery)	16kg/35lbs
Type of Use	Pavement and indoor
Braking System	Electro-magnetic
Airline friendly	Yes

\*Specifications may vary due to manufacturing tolerances



**7 Year Battery + 2 Year Scooter Warranty\*\***

\*\*terms and conditions apply



Scan for more  
information

eFOLDi 

**For more information...**

Call us today on  
**0800 999 8123**

Or email us at:  
**[sales@efoldi.com](mailto:sales@efoldi.com)**

Your local eFOLDi retailer:



SunTech UK Limited t/a eFOLDi, 25 Ormside Way,  
Holmethorpe Industrial Estate, Redhill, Surrey RH1 2LW.