



7 Year
Battery + 2
Year Chair
Warranty**

Enjoy hassle-free travel again and make every journey effortless

The eFOLDi Navigator is one of the world's lightest, most compact folding powered wheelchairs.

Weighing just 12.4kg*, the ultra-portable Navigator folds and unfolds in just a few seconds so is easy to store and transport, making it the perfect travel companion.

No complex assembly is required – just unfold and go.



Lightweight
at only
12.4kg/27lbs*




Top speed
4mph/6kph



Range
up to
8miles/13km



Warranty
Backed by
full support &
warranty


Ideal For: Those with mobility issues who prioritise extra comfort and ease, looking for a power chair that supports their everyday travel needs with added convenience.



The lighter way to empower your life

The perfect blend of uncompromising style, portability and versatility

Integrated attendant control
switch easily between self and assisted driving

Simple, ergonomic joystick-style controls
for smooth, intuitive handling

Automotive grade Lithium-ion battery
with 8 miles/13km range, removable for on/off-board charging

Durable frame
made from tough magnesium alloy to handle everyday use with ease

Anti-tip wheels
for enhanced safety and reassurance

Solid tyres
with no pumping requirement and multi-terrain suitability

Powerful low-vibration motor
with whisper-quiet performance

Cushion
Removable, breathable seat cushion for optimum comfort

Electro-magnetic brake system
for maximum safety/control

Adjustable footrest
for comfort and accessibility across a wide range of users

Compact
Folds into a compact size within seconds, so easy to store and transport

Ultra-lightweight
at only 12.4kg/27lbs*



"Being so lightweight and portable, means my eFOLDi really helps me and my wife make the most of our holiday getaways, at home and abroad."

James, Verified Buyer

Lightweight, comfortable, quick and simple to fold... the Navigator combines it all.

- Easy to use, ultra-portable folding powered wheelchair
- Weighing only 12.4kg*, it's one of the world's lightest, folding powerchairs available
- Folds in seconds, fits into most car boots and is 'Airline-safe'
- Ergonomically designed user controls, plus integral attendant control for total peace of mind
- Advanced, high performance battery technology with on- or off-board charging options

Accessories

Make every journey easier with thoughtful extras designed for your comfort and convenience – from additional cushions and attachments to spare batteries and chargers.



* without battery/accessories
**terms and conditions apply



COMPASS

LITE

EXPLORER

NAVIGATOR

The eFOLDi Navigator offers superior levels of portability, convenience and versatility

As a recognised leader in the design and manufacture of portable scooters and powered wheelchairs, eFOLDi strives to help people rediscover their freedom and independence.

eFOLDi's award-winning mobility solutions are not only easy to store and transport, they help make travel simple and achievable – whether that's by air, land or sea.

Why eFOLDi stands apart




Lightweight
Carry and lift with ease



Warranty
Backed by full support & warranty



Compact
Folds down to 'hold' luggage size



Robust
Peace of mind with strong, durable materials

Innovation

Our products are manufactured using the highest quality materials to ensure they're robust enough to withstand the toughest of conditions and terrain. Our frames are made from high-grade aluminium or magnesium alloys, our high-specification brushless motors deliver exceptional performance, and our batteries use advanced lithium-ion technology, so they're not only durable, they offer outstanding levels of reliability.

Navigator Specifications*

Maximum Speed	4mph / 6kph
Maximum Range	8 miles / 13km
Maximum User Weight	19 Stone / 120kg / 265lbs
Climbing Ability	1:10 Gradient
Turning Radius	120cm (40°)
Folded Dimensions	34 x 59 x 89cm 13 x 23 x 35 inches
Unfolded Dimensions	108 x 59 x 89cm 43 x 23 x 35 inches
Weight (excl. battery)	12.4kg/27lbs
Weight (incl. battery)	14kg/31lbs
Type of Use	Pavement and indoor
Braking System	Electro-magnetic
Airline friendly	Yes

*Specifications may vary due to manufacturing tolerances



7 Year Battery + 2 Year Chair Warranty**

**terms and conditions apply



Scan for more information

eFOLDi

For more information...

Call us today on
0800 999 8123

Or email us at:
sales@efoldi.com



SunTech UK Limited t/a eFOLDi, 25 Ormside Way,
Holmetherpe Industrial Estate, Redhill, Surrey RH1 2LW.

eFOLDi

eFOLDi Navigator

Lightweight mobility solutions

For wherever life takes you

www.efoldi.com



Scan for more information