

# Lightweight mobility solutions

For wherever life takes you

[www.efoldi.com](http://www.efoldi.com)



Scan for more  
information



**7 Year  
Battery + 2  
Year Chair  
Warranty\*\***

## Enjoy hassle-free travel again and make every journey effortless

The eFOLDi Navigator is one of the world's lightest, most compact folding powered wheelchairs.

Weighing just 12.4kg\*, the ultra-portable Navigator folds and unfolds in just a few seconds so is easy to store and transport, making it the perfect travel companion.

No complex assembly is required – just unfold and go.



**Lightweight**  
at only  
12.4kg/27lbs\*



**Top speed**  
4mph/6kph



**Range**  
up to  
8miles/13km



**Warranty**  
Backed by  
full support &  
warranty



**Ideal For:** Those with mobility issues who prioritise extra comfort and ease, looking for a power chair that supports their everyday travel needs with added convenience.




\* without battery/accessories

\*\*terms and conditions apply

# The lighter way to empower your life

The perfect blend of uncompromising style, portability and versatility





"Being so lightweight and portable, means my eFOLDi really helps me and my wife make the most of our holiday getaways, at home and abroad."

James, Verified Buyer 

## Lightweight, comfortable, quick and simple to fold... the Navigator combines it all.

- Easy to use, ultra-portable folding powered wheelchair
- Weighing only 12.4kg\*, it's one of the world's lightest, folding powerchairs available
- Folds in seconds, fits into most car boots and is 'Airline-safe'
- Ergonomically designed user controls, plus integral attendant control for total peace of mind
- Advanced, high performance battery technology with on- or off-board charging options

## Accessories

Make every journey easier with thoughtful extras designed for your comfort and convenience – from additional cushions and attachments to spare batteries and chargers.



## The eFOLDi Navigator offers superior levels of portability, convenience and versatility

As a recognised leader in the design and manufacture of portable scooters and powered wheelchairs, eFOLDi strives to help people rediscover their freedom and independence.

eFOLDi's award-winning mobility solutions are not only easy to store and transport, they help make travel simple and achievable – whether that's by air, land or sea.

COMPASS



LITE



EXPLORER



NAVIGATOR



## Why eFOLDi stands apart



### Lightweight

Carry and lift with ease




### Warranty

Backed by full support & warranty



### Compact

Folds down to 'hold' luggage size



### Robust

Peace of mind with strong, durable materials



## Innovation

Our products are manufactured using the highest quality materials to ensure they're robust enough to withstand the toughest of conditions and terrain. Our frames are made from high-grade aluminium or magnesium alloys, our high-specification brushless motors deliver exceptional performance, and our batteries use advanced lithium-ion technology, so they're not only durable, they offer outstanding levels of reliability.

## Navigator Specifications\*

Maximum Speed	4mph / 6kph
Maximum Range	8 miles / 13km
Maximum User Weight	19 Stone / 120kg / 265lbs
Climbing Ability	1:10 Gradient
Turning Radius	120cm (40*)
Folded Dimensions	34 x 59 x 89cm 13 x 23 x 35 inches
Unfolded Dimensions	108 x 59 x 89cm 43 x 23 x 35 inches
Weight (excl. battery)	12.4kg/27lbs
Weight (incl. battery)	14kg/31lbs
Type of Use	Pavement and indoor
Braking System	Electro-magnetic
Airline friendly	Yes

\*Specifications may vary due to manufacturing tolerances



**7 Year Battery + 2 Year Chair Warranty\*\***

\*\*terms and conditions apply



Scan for more  
information

eFOLDi 

**For more information...**

Call us today on **0800 999 8123** Or email us at: **[sales@efoldi.com](mailto:sales@efoldi.com)**

Your local eFOLDi retailer:



SunTech UK Limited t/a eFOLDi, 25 Ormside Way,  
Holmethorpe Industrial Estate, Redhill, Surrey RH1 2LW.