



Lightweight mobility solutions

For wherever life takes you

www.efoldi.com



Wherever life takes you, we'll be right here, making sure you have the mobility and freedom to enjoy every moment.

My name is Sumi, and welcome to eFOLDi – Home to some of the world's lightest and most compact folding mobility scooters and power chair.

We strive to help people rediscover their freedom and independence. When my father lost his mobility in a work-place accident some years ago, he simply couldn't find a mobility scooter anywhere, that was easy to use, safe and comfortable – but also lightweight and portable. Undeterred, my father set about designing his own scooter – and that is how the first eFOLDi mobility scooter came to be born! Since then, eFOLDi has developed into a thriving family business, designing and manufacturing innovative mobility products that help make travel simple, achievable and most importantly, fun!

The eFOLDi scooter was also a proud winner of the prestigious 'British Invention Of The Year' award in 2016. At the Virgin Media VOOM competition in 2016, legendary innovator and entrepreneur, Sir Richard Branson, even named the eFOLDi scooter "folding magic".

Join us in making travel easy and enjoyable again, and let an eFOLDi be your perfect travel companion.

Sumi Wang
Founder

COMPASS



CLIQ



LITE



EXPLORER



NAVIGATOR



Why eFOLDi stands apart



Lightweight

Carry and lift
with ease




Warranty

Backed by full
support & warranty



Compact

Folds down to 'hold'
luggage size



Robust

Peace of mind with
strong, durable
materials



"I like the independence the eFOLDi gives me, I'm not having to wait for other people - I can go where I want when I want."

Maria, Verified Buyer 



eFOLDi Compass

Stable, Comfortable & Ready for Every Journey

Our first 4-wheel, bi-folding scooter, the new Compass is our easiest, fastest folding scooter yet—designed to make your journeys easier and more enjoyable. Weighing just 11kg (without the battery), the Compass is lightweight, stable and has a surprisingly tight turning circle for a 4-wheel scooter, making it nimble to get around on. Its compact design allows for effortless travel, ensuring you can continue exploring, socialising, and enjoying life on your terms.



Lightweight

at only
11kg/24lbs*



Top speed

4mph/6kph



Range

10miles/16km



Warranty

Backed by
full support &
warranty

*without battery/accessories



Travel in comfort

Optional armrests
for extra comfort



Compact size

Folds with ease



Lightweight

at only 11kg/24lbs*



Robust

Solid,
low-maintenance
rear tyres



Secure-Turn™ technology

helps enhance stability
and control by detecting
changes in movement and
orientation



Airline friendly battery

10 miles/16km per charge.
Allows on-or-off-board
charging



**7 Year
Battery + 2
Year Scooter**
Warranty****

"It's so clever —
the scooter adjusts
automatically when
I turn, which makes
driving feel safe and
smooth."

Ros, Verified Buyer



*without battery/accessories



eFOLDi Cliq


Comfort and portability, clicked together.

The eFOLDi Cliq has four parts each under 10 kg, designed with EZ-Cliq connections for quick, independent assembly. Enjoy longer trips with a smooth, supportive ride. For those who value stability, a comfortable seat, and a planted wheelbase, with the flexibility to take the scooter apart for easy transport.



Lightweight

31kg total weight
with 4 easy-lift
parts



Robust Build

Carries
users up to
120kg/265lbs/
19 Stone



Range

Up to
9miles/14km



Warranty

Backed by
full support
& warranty



**7 Year
Battery + 2
Year Scooter
Warranty****

"The comfortable seat means I can enjoy longer days out, and because it swivels, getting on and off is so much easier."

Les, Verified Buyer



Travel in comfort

Flip-up armrests and practical storage basket for real-world convenience



Lightweight

31kg total weight with 4 easy-lift parts



Personalisation

Interchangeable coloured accents as standard to customise your ride



Robust

Carries users up to 120kg/265lbs/19 Stone



Compact size

individual parts weigh <10kg with easy-Click connections for effortless assembly and transport



eFoldi Yellow



Moonlight Silver



Harbour Blue





eFOLDi Lite

Light, Compact & Ready for Adventure

Our lightest most compact three-wheeled scooter. Weighing just 14kg (without the battery), designed for ultimate portability, fitting easily into a small car or on public transport, so you can take it anywhere with ease.



Lightweight

at only
14kg/31lbs*



Top speed

4mph/6kph



Range

10miles/16km



Warranty

Backed by
full support &
warranty

*without battery/accessories



Travel in comfort

Optional armrests
for extra comfort



Compact size

Folds with ease



Lightweight

at only 14kg/31lbs*



Robust

Solid,
low-maintenance
rear tyres



Speed reducer

Automatically slows
down when turning



Airline friendly battery

10 miles/16km per charge.
Allows on-or-off-board
charging

**7 Year
Battery + 2
Year Scooter
Warranty****



"I can go out when I want
and don't have to worry
if there is some walking
involved."

Diane, Verified Buyer

*without battery/accessories **for purchases made after 01/01/25



eFOLDi Explorer



Robust, Versatile & Ready for Any Terrain

Our multi award-winning, road-legal* mobility scooter that combines lightweight design with powerful performance. Travel up to 8mph and enjoy a range of up to 14 miles, all while gaining a new sense of freedom.



Lightweight
at only
17kg/38lbs**



Top speed
8mph/12kph



Range
up to
14miles/22km



Warranty
Backed by
full support &
warranty

*Approved for UK roads only (Hazard light must be used on dual carriageways)

**without battery/accessories



**7 Year
Battery + 2
Year Scooter
Warranty****

"With the eFOLDi I'm able to visit places I haven't been able to get to for a while. I'm very happy to tell people – buy an eFOLDi."

Peter, Verified Buyer ✓



Compact size
folds with ease



Fully UK road legal
8mph/12kph
top speed



Lightweight
at only 17kg/38lbs*



**Airline friendly
battery**
14 miles/22km per charge. Allows
on-or-off-board charging



Easy to use
controls



Travel in comfort
suspension/
pneumatic tyres

*without battery/accessories **for purchases made after 01/01/25

eFOLDi Navigator

Quiet, Lightweight & Effortless Mobility

Our ultra-light, whisper-quiet powerchair, designed for comfort and effortless everyday travel. Weighing just 12kg (without the battery), it folds in seconds to store or transport and with a discreet attendant control as standard, it offers ease and convenience for both users and their companions.



Lightweight
at only
12kg/26lbs*

*without battery/accessories



Top speed
4mph/6kph



Range
8miles/13km



Warranty
Backed by
full support &
warranty





Compact size
folds with ease



Travel in comfort
durable,
breathable fabric



Handy controls
for attendant to drive



**Airline friendly
battery**
8 miles/13km per charge. Allows
on-or-off-board charging



Easy to use
joystick and control pad

**7 Year Battery
+ 2 Year Chair
Warranty****



Lightweight
at only 12kg/26lbs*



"With an eFOLDi you can just fold it up, and take it with you on buses, trains, or planes. Perfect for me to get out and about with my family and even go on holiday abroad!"

Anna, Verified Buyer  13

*without battery/accessories **for purchases made after 01/01/25



	eFOLDi Compass	eFOLDi Cliq
Maximum Speed	4mph / 6kph	4mph / 6kph
Maximum Range	10 Miles / 16km	9 Miles / 15km
Maximum User Weight	19 Stone / 120kg / 265lbs	19 Stone / 120kg / 265lbs
Climbing Ability	1:10 Gradient	1:10 Gradient
Turning Radius	135cm (53*)	121cm (48*)
Folded Dimensions (L/W/H)	36 x 17 x 16 inches 92 x 43 x 40cm	N/A
Unfolded Dimensions (L/W/H)	36 x 17 x 38 inches 92 x 43 x 96cm	39 x 24 x 34.6 inches 100 x 61 x 88cm
Individual Part Weight	N/A	3.6kg/7.9lbs 9.3kg/20.5lbs 9.8kg/21.6lbs 7.3kg/15.9lbs
Weight (incl. battery)	13kg / 29lbs	31kg / 69lbs
Weight (excl. battery)	11kg / 25lbs	28kg / 61lbs
Type of Use	Pavement and Indoor	Pavement and Indoor
Braking System	Electro-Magnetic	Electro-Magnetic
Airline friendly	Yes	Yes

eFOLDi Lite	eFOLDi Explorer	eFOLDi Navigator
4mph / 6kph	8mph / 12kph	4mph / 6kph
10 miles / 16km	14 miles / 22km	8 miles / 13km
19 Stone / 120kg / 265lbs	19 Stone / 120kg / 265lbs	21 Stone / 135kg / 297lbs
1:10 Gradient	1:6 Gradient	1:10 Gradient
100cm (40*)	100cm (40*)	120cm (40*)
27.5 x 22 x 15 inches 70 x 55 x 37cm	26 x 22 x 16 inches 66 x 56 x 40cm	13 x 23 x 35 inches 34 x 59 x 89cm
46 x 22 x 37 inches 117 x 55 x 94cm	48 x 22 x 33 inches 121 x 56 x 83cm	43 x 23 x 35 inches 108 x 59 x 89cm
N/A	N/A	N/A
16kg/35lbs	19kg/42lbs	14kg/31lbs
14kg/31lbs	17kg/38lbs	12.4kg/27lbs
Pavement and Indoor	Road (UK legal)/Pavement, Indoor	Pavement and Indoor
Electro-Magnetic	Mechanical Hand Brake	Electro-Magnetic
Yes	Yes	Yes

Specifications may vary due to manufacturing tolerances.



// I've always been active, especially after my career as a firefighter, so suddenly losing that freedom was frustrating. Getting my eFOLDi Explorer changed everything, it gave me my freedom back. It lives in the car boot, so I can make the most of different locations and join in with walks, even when we're not near a pavement. I'm not embarrassed to use it at all; people are always amazed when they see how it folds up."

Denis, Verified Buyer 

Accessories

Scooter Lifter

Our lightweight scooters are designed for travel, but for those with limited strength or restricted movement, even a few kilos can feel like a challenge. Our electronic lifter comes is compact, portable, and powered directly from your scooter battery, it gives you extra loading support when you need it.



*for scooters only

All Product Accessories



Batteries



High Vis Bag



Hi-Vis Vest



Chargers

Lite



Armrest



Red Rear Seat Pad



Storage Case



Soft Grip Handle Extension

Red Main Seat Pad



Comfort Cushion



Protection Cover



Basket




Footrest

Navigator



Protection Cover



Extension Bar



Cushion



Under Seat Bag



Cup Holder

Explorer



Knobbly Tyre




Tyre Inflator



Protection Cover



Comfort Cushion



Soft Grip Handle Extension




Foot Peg Ext (Single)

Compass



Storage Case



Protection Cover



Steadytrred Wheels



Soft Grip Handle Extension



Basket



Armrests



Innovation

Recognised as leaders in the design and manufacture of ultra-portable scooters and power chairs, we strive to help people rediscover their freedom and independence.

Our products are manufactured using the highest quality materials to ensure they're robust enough to withstand the toughest of conditions and terrain. Our scooter frames are made from high-grade aluminium or magnesium alloys, our high-specification brushless motors deliver exceptional performance, and our batteries use advanced lithium-ion technology, so they're not only durable, they offer outstanding levels of reliability.



Support

All eFOLDi products come with a 7-year battery warranty and 2-year scooter/chair warranty. In the unlikely event of a problem, our UK support centre and servicing team are always on hand to help. They are trained to diagnose and rectify any problems, including warranty repairs and any ongoing service and maintenance requirements.

Home Demos

At eFOLDi, we believe mobility should be simple, stress-free, and tailored to you. That's why we bring the showroom experience straight to your doorstep! No pressure, no hassle—just friendly service and a chance to try out our scooters and power chairs in your own space, at a time that works for you.



eFOLDi 



“For years I sat inside...
and then came this
life changer!”

Tracey, Verified Buyer 



EXPLORER



Scan for more
information

Your local eFOLDi retailer:



SunTech UK Limited t/a eFOLDi, 25 Ormside Way,
Holmethorpe Industrial Estate, Redhill, Surrey RH1 2LW.

For more information...

Call us today on
0800 999 8123

Or email us at:
sales@efoldi.com