



# Lightweight mobility solutions

For wherever life takes you

[www.efoldi.com](http://www.efoldi.com)



# Wherever life takes you, we'll be right here, making sure you have the mobility and freedom to enjoy every moment.

**My name is Sumi, and welcome to eFOLDi – Home to some of the world's lightest and most compact folding mobility scooters and power chair.**

We strive to help people rediscover their freedom and independence. When my father lost his mobility in a work-place accident some years ago, he simply couldn't find a mobility scooter anywhere, that was easy to use, safe and comfortable – but also lightweight and portable. Undeterred, my father set about designing his own scooter – and that is how the first eFOLDi mobility scooter came to be born! Since then, eFOLDi has developed into a thriving family business, designing and manufacturing innovative mobility products that help make travel simple, achievable and most importantly, fun!

The eFOLDi scooter was also a proud winner of the prestigious 'British Invention Of The Year' award in 2016. At the Virgin Media VOOM competition in 2016, legendary innovator and entrepreneur, Sir Richard Branson, even named the eFOLDi scooter "folding magic".

Join us in making travel easy and enjoyable again, and let an eFOLDi be your perfect travel companion.

*Sumi Wang*  
**Founder**

COMPASS



CLIQ



LITE



EXPLORER



NAVIGATOR



# Why eFOLDi stands apart



## Lightweight

Carry and lift  
with ease




## Warranty

Backed by full  
support & warranty



## Compact

Folds down to 'hold'  
luggage size



## Robust

Peace of mind with  
strong, durable  
materials



"I like the independence the eFOLDi gives me, I'm not having to wait for other people - I can go where I want when I want."

Maria, Verified Buyer 



## eFOLDi Compass

### Stable, Comfortable & Ready for Every Journey

Our first 4-wheel, bi-folding scooter, the new Compass is our easiest, fastest folding scooter yet—designed to make your journeys easier and more enjoyable. Weighing just 11kg (without the battery), the Compass is lightweight, stable and has a surprisingly tight turning circle for a 4-wheel scooter, making it nimble to get around on. Its compact design allows for effortless travel, ensuring you can continue exploring, socialising, and enjoying life on your terms.



#### Lightweight

at only  
11kg/24lbs\*



#### Top speed

4mph/6kph



#### Range

10miles/16km



#### Warranty

Backed by  
full support &  
warranty

\*without battery/accessories





### Travel in comfort

Optional armrests  
for extra comfort



### Compact size

Folds with ease



### Lightweight

at only 11kg/24lbs\*



### Robust

Solid,  
low-maintenance  
rear tyres



### Secure-Turn™ technology

helps enhance stability  
and control by detecting  
changes in movement and  
orientation



### Airline friendly battery

10 miles/16km per charge.  
Allows on-or-off-board  
charging



**7 Year  
Battery + 2  
Year Scooter\*\*  
Warranty\*\***

"It's so clever —  
the scooter adjusts  
automatically when  
I turn, which makes  
driving feel safe and  
smooth."

**Ros**, Verified Buyer



\*without battery/accessories



## eFOLDi Cliq


### Comfort and portability, clicked together.

The eFOLDi Cliq has four parts each under 10 kg, designed with EZ-Cliq connections for quick, independent assembly. Enjoy longer trips with a smooth, supportive ride. For those who value stability, a comfortable seat, and a planted wheelbase, with the flexibility to take the scooter apart for easy transport.



#### Lightweight

31kg total weight  
with 4 easy-lift  
parts



#### Robust Build

Carries  
users up to  
120kg/265lbs/  
19 Stone



#### Range

Up to  
9miles/14km



#### Warranty

Backed by  
full support  
& warranty



**7 Year  
Battery + 2  
Year Scooter  
Warranty\*\***

"The comfortable seat means I can enjoy longer days out, and because it swivels, getting on and off is so much easier."

Les, Verified Buyer



### Travel in comfort

Flip-up armrests and practical storage basket for real-world convenience



### Lightweight

31kg total weight with 4 easy-lift parts



### Personalisation

Interchangeable coloured accents as standard to customise your ride



### Robust

Carries users up to 120kg/265lbs/19 Stone



### Compact size

individual parts weigh <10kg with easy-Click connections for effortless assembly and transport



eFoldi Yellow



Moonlight Silver



Harbour Blue





## eFOLDi Lite

# Light, Compact & Ready for Adventure

Our lightest most compact three-wheeled scooter. Weighing just 14kg (without the battery), designed for ultimate portability, fitting easily into a small car or on public transport, so you can take it anywhere with ease.



### Lightweight

at only  
14kg/31lbs\*



### Top speed

4mph/6kph



### Range

10miles/16km



### Warranty

Backed by  
full support &  
warranty

\*without battery/accessories





### Travel in comfort

Optional armrests  
for extra comfort



### Compact size

Folds with ease



### Lightweight

at only 14kg/31lbs\*



### Robust

Solid,  
low-maintenance  
rear tyres



### Speed reducer

Automatically slows  
down when turning



### Airline friendly battery

10 miles/16km per charge.  
Allows on-or-off-board  
charging

**7 Year  
Battery + 2  
Year Scooter  
Warranty\*\***



"I can go out when I want  
and don't have to worry  
if there is some walking  
involved."

**Diane**, Verified Buyer

\*without battery/accessories \*\*for purchases made after 01/01/25



## eFOLDi Explorer



### Robust, Versatile & Ready for Any Terrain

Our multi award-winning, road-legal\* mobility scooter that combines lightweight design with powerful performance. Travel up to 8mph and enjoy a range of up to 14 miles, all while gaining a new sense of freedom.



**Lightweight**  
at only  
17kg/38lbs\*\*



**Top speed**  
8mph/12kph



**Range**  
up to  
14miles/22km



**Warranty**  
Backed by  
full support &  
warranty

\*Approved for UK roads only (Hazard light must be used on dual carriageways)

\*\*without battery/accessories



**7 Year  
Battery + 2  
Year Scooter  
Warranty\*\***

"With the eFOLDi I'm able to visit places I haven't been able to get to for a while. I'm very happy to tell people – buy an eFOLDi."

**Peter**, Verified Buyer ✓



**Compact size**  
folds with ease



**Fully UK road legal**  
8mph/12kph  
top speed



**Lightweight**  
at only 17kg/38lbs\*



**Airline friendly  
battery**  
14 miles/22km per charge. Allows  
on-or-off-board charging



**Easy to use**  
controls



**Travel in comfort**  
suspension/  
pneumatic tyres

\*without battery/accessories \*\*for purchases made after 01/01/25

# eFOLDi Navigator

## Quiet, Lightweight & Effortless Mobility

Our ultra-light, whisper-quiet powerchair, designed for comfort and effortless everyday travel. Weighing just 12kg (without the battery), it folds in seconds to store or transport and with a discreet attendant control as standard, it offers ease and convenience for both users and their companions.



**Lightweight**  
at only  
12kg/26lbs\*

\*without battery/accessories



**Top speed**  
4mph/6kph



**Range**  
8miles/13km



**Warranty**  
Backed by  
full support &  
warranty







**Compact size**  
folds with ease



**Travel in comfort**  
durable,  
breathable fabric



**Handy controls**  
for attendant to drive



**Airline friendly  
battery**  
8 miles/13km per charge. Allows  
on-or-off-board charging



**Easy to use**  
joystick and control pad

**7 Year Battery  
+ 2 Year Chair  
Warranty\*\***



**Lightweight**  
at only 12kg/26lbs\*



"With an eFOLDi you can just fold it up, and take it with you on buses, trains, or planes. Perfect for me to get out and about with my family and even go on holiday abroad!"

**Anna, Verified Buyer**  13

\*without battery/accessories \*\*for purchases made after 01/01/25



|                             | eFOLDi Compass                        | eFOLDi Cliq                                                   |
|-----------------------------|---------------------------------------|---------------------------------------------------------------|
| Maximum Speed               | 4mph / 6kph                           | 4mph / 6kph                                                   |
| Maximum Range               | 10 Miles / 16km                       | 9 Miles / 15km                                                |
| Maximum User Weight         | 19 Stone /<br>120kg / 265lbs          | 19 Stone /<br>120kg / 265lbs                                  |
| Climbing Ability            | 1:10 Gradient                         | 1:10 Gradient                                                 |
| Turning Radius              | 135cm (53*)                           | 121cm (48*)                                                   |
| Folded Dimensions (L/W/H)   | 36 x 17 x 16 inches<br>92 x 43 x 40cm | N/A                                                           |
| Unfolded Dimensions (L/W/H) | 36 x 17 x 38 inches<br>92 x 43 x 96cm | 39 x 24 x 34.6 inches<br>10 x 61 x 88cm                       |
| Individual Part Weight      | N/A                                   | 3.6kg/7.9lbs   9.3kg/20.5lbs<br>9.8kg/21.6lbs   7.3kg/15.9lbs |
| Weight (incl. battery)      | 13kg / 29lbs                          | 28kg / 61lbs                                                  |
| Weight (excl. battery)      | 11kg / 25lbs                          | 31kg / 69lbs                                                  |
| Type of Use                 | Pavement and Indoor                   | Pavement and Indoor                                           |
| Braking System              | Electro-Magnetic                      | Electro-Magnetic                                              |
| Airline friendly            | Yes                                   | Yes                                                           |

| eFOLDi Lite                             | eFOLDi Explorer                        | eFOLDi Navigator                       |
|-----------------------------------------|----------------------------------------|----------------------------------------|
| 4mph / 6kph                             | 8mph / 12kph                           | 4mph / 6kph                            |
| 10 miles / 16km                         | 14 miles / 22km                        | 8 miles / 13km                         |
| 19 Stone /<br>120kg / 265lbs            | 19 Stone /<br>120kg / 265lbs           | 21 Stone /<br>135kg / 297lbs           |
| 1:10 Gradient                           | 1:6 Gradient                           | 1:10 Gradient                          |
| 100cm (40*)                             | 100cm (40*)                            | 120cm (40*)                            |
| 27.5 x 22 x 15 inches<br>70 x 55 x 37cm | 26 x 22 x 16 inches<br>66 x 56 x 40cm  | 13 x 23 x 35 inches<br>34 x 59 x 89cm  |
| 46 x 22 x 37 inches<br>117 x 55 x 94cm  | 48 x 22 x 33 inches<br>121 x 56 x 83cm | 43 x 23 x 35 inches<br>108 x 59 x 89cm |
| N/A                                     | N/A                                    | N/A                                    |
| 16kg/35lbs                              | 19kg/42lbs                             | 14kg/31lbs                             |
| 14kg/31lbs                              | 17kg/38lbs                             | 12.4kg/27lbs                           |
| Pavement and Indoor                     | Road (UK legal)/Pavement,<br>Indoor    | Pavement and Indoor                    |
| Electro-Magnetic                        | Mechanical Hand Brake                  | Electro-Magnetic                       |
| Yes                                     | Yes                                    | Yes                                    |

Specifications may vary due to manufacturing tolerances.



I've always been active, especially after my career as a firefighter, so suddenly losing that freedom was frustrating. Getting my eFOLDi Explorer changed everything, it gave me my freedom back. It lives in the car boot, so I can make the most of different locations and join in with walks, even when we're not near a pavement. I'm not embarrassed to use it at all; people are always amazed when they see how it folds up."

Denis, Verified Buyer 

# Accessories

## Scooter Lifter

Our lightweight scooters are designed for travel, but for those with limited strength or restricted movement, even a few kilos can feel like a challenge. Our electronic lifter comes is compact, portable, and powered directly from your scooter battery, it gives you extra loading support when you need it.



\*for scooters only

## All Product Accessories



Batteries

High Vis Bag

Chargers

Hi-Vis Vest

## Lite



Armrest

Red Rear Seat Pad

Storage Case

Red Main Seat Pad

Soft Grip Handle Extension

Comfort Cushion

Protection Cover

Footrest


Basket



## Navigator



Protection Cover



Extension Bar



Cushion



Under Seat Bag



Cup Holder

## Explorer



Knobbly Tyre




Tyre Inflator



Protection Cover



Comfort Cushion



Soft Grip Handle Extension




Foot Peg Ext (Single)

## Compass



Storage Case



Protection Cover



Steadytred Wheels



Soft Grip Handle Extension



Basket



Armrests



## Innovation

**Recognised as leaders in the design and manufacture of ultra-portable scooters and power chairs, we strive to help people rediscover their freedom and independence.**

Our products are manufactured using the highest quality materials to ensure they're robust enough to withstand the toughest of conditions and terrain. Our scooter frames are made from high-grade aluminium or magnesium alloys, our high-specification brushless motors deliver exceptional performance, and our batteries use advanced lithium-ion technology, so they're not only durable, they offer outstanding levels of reliability.



## Support

All eFOLDi products come with a 7-year battery warranty and 2-year scooter/chair warranty. In the unlikely event of a problem, our UK support centre and servicing team are always on hand to help. They are trained to diagnose and rectify any problems, including warranty repairs and any ongoing service and maintenance requirements.

## Home Demos

At eFOLDi, we believe mobility should be simple, stress-free, and tailored to you. That's why we bring the showroom experience straight to your doorstep! No pressure, no hassle—just friendly service and a chance to try out our scooters and power chairs in your own space, at a time that works for you.



eFOLDi 



“For years I sat inside...  
and then came this  
life changer!”

Tracey, Verified Buyer 



**EXPLORER**



Scan for more  
information

Your local eFOLDi retailer:



SunTech UK Ltd, Unit H, The Courtyard Business Centre,  
Dovers Farm, Lonesome lane, Reigate, Surrey RH2 7QT.

**For more information...**

Call us today on  
**0800 999 8123**

Or email us at:  
**[sales@efoldi.com](mailto:sales@efoldi.com)**