



# eFOLDi<sup>®</sup> Information for Airlines and Handlers

## Device information

|                  |                              |
|------------------|------------------------------|
| Manufacturer:    | eFOLDi                       |
| Model:           | Navigator                    |
| Serial Number:   |                              |
| Device Value:    | £2,699                       |
| Maximum Weight:  | 14kg                         |
| Weight Unloaded: | 12.4kg                       |
| Height:          | 34cm folded / 108cm unfolded |
| Width:           | 59cm folded / 59cm unfolded  |
| Length:          | 89cm folded / 89cm unfolded  |

## Manoeuvrability & adjustment features

- Driving with controller
- Push in freewheel only
- Back folds down
- Back reclines
- Seat tilts
- Release lever marked

## Battery Information

|                 |             |
|-----------------|-------------|
| Battery type:   | Lithium-ion |
| Quantity:       | one         |
| Watt-hour (Wh): | 264.6 Wh    |



# eFOLDi<sup>®</sup> Information for Airlines and Handlers

## Removeable parts

| Component         | Remain on device | Cabin |
|-------------------|------------------|-------|
| Seat Cushion      | ✓                |       |
| Head Support      | N/A              |       |
| Arm Supports      | ✓                |       |
| Leg Supports      | N/A              |       |
| Control Device    | ✓                |       |
| Secondary Control | ✓                |       |
| Backrest          | ✓                |       |
| Tray              | N/A              |       |
| Belt/Straps       | ✓                |       |
| Other             | N/A              |       |

## Passenger Contact

Name: \_\_\_\_\_

Email: \_\_\_\_\_

Mobile: \_\_\_\_\_

Flight Number: \_\_\_\_\_

**CRITICAL: Lift only from base/frame - NEVER from armrests**