



eFOLDi[®] Information for Airlines and Handlers

Device information

Manufacturer:	eFOLDi
Model:	Compass
Serial Number:	
Device Value:	£2,499
Maximum Weight:	13kg
Weight Unloaded:	11kg
Height:	40cm folded / 96cm unfolded
Width:	43cm folded / 43cm unfolded
Length:	92cm folded / 92cm unfolded

Manoeuvrability & adjustment features

- Driving with controller
- Push in freewheel only
- Back folds down
- Back reclines
- Seat tilts
- Release lever marked

Battery Information

Battery type:	Lithium-ion
Quantity:	one
Watt-hour (Wh):	264.6 Wh



eFOLDi[®] Information for Airlines and Handlers

Removeable parts

Component	Remain on device	Cabin
Seat Cushion	✓	
Head Support	N/A	
Arm Supports	✓	
Leg Supports	N/A	
Control Device	✓	
Secondary Control	N/A	
Backrest	N/A	
Tray	N/A	
Belt/Straps	N/A	
Other	N/A	

Passenger Contact

Name: _____

Email: _____

Mobile: _____

Flight Number: _____

CRITICAL: Lift only from base/frame - NEVER from armrests