



# eFOLDi<sup>®</sup> Information for Airlines and Handlers

## Device information

Manufacturer:	eFOLDi
Model:	Lite
Serial Number:	
Device Value:	£2,699
Maximum Weight:	16kg
Weight Unloaded:	14kg
Height:	70cm folded / 117cm unfolded
Width:	55cm folded / 55cm unfolded
Length:	94cm folded / 94cm unfolded

## Manoeuvrability & adjustment features

- Driving with controller
- Push in freewheel only
- Back folds down
- Back reclines
- Seat tilts
- Release lever marked

## Battery Information

Battery type:	Lithium-ion
Quantity:	one
Watt-hour (Wh):	264.6 Wh



# eFOLDi<sup>®</sup> Information for Airlines and Handlers

## Removeable parts

Component	Remain on device	Cabin
Seat Cushion	✓	
Head Support	N/A	
Arm Supports	✓	
Leg Supports	N/A	
Control Device	✓	
Secondary Control	N/A	
Backrest	N/A	
Tray	N/A	
Belt/Straps	N/A	
Other	N/A	

## Passenger Contact

**Name:** \_\_\_\_\_

**Email:** \_\_\_\_\_

**Mobile:** \_\_\_\_\_

**Flight Number:** \_\_\_\_\_

**CRITICAL: Lift only from base/frame - NEVER from armrests**