



7 Year Battery + 2 Year Scooter Warranty**

Comfort and portability, clicked together.

The eFOLDi Cliq has four parts each under 10 kg, designed with EZ-Cliq connections for quick, independent assembly.

Enjoy longer trips with a smooth, supportive ride.

For those who value stability, a comfortable seat, and a planted wheelbase, with the flexibility to take the scooter apart for easy transport.



Lightweight

31kg total weight with 4 easy-lift parts



Robust Build

Carries users up to 120kg/265lbs/19 Stone



Range

up to 9miles/14km



Warranty

Backed by full support & warranty*

*terms and conditions apply

Ideal For: individuals with moderate mobility issues who need extra support and stability, seeking a portable scooter that offers increased confidence, security, and ease of use for getting around.



Your Confident Companion on Every Journey

Designed for comfort, stability, and everyday independence.



“With the eFOLDi I’m able to visit places I haven’t been able to get to for a while. I’m very happy to tell people – buy an eFOLDi.”

Peter, Verified Buyer

Clicks Together

- Modular design: easy assembly and transportation
- Lightweight parts: <10kg with Easy-Cliq connections
- Snap-together simplicity: hassle-free setup
- Fits everywhere: car and home storage without stress

Clicks With You

- Comfort-first seating: adjustable, supportive seat with all-round suspension
- Confident Controls: WigWag levers, ergonomic delta handlebars, simple speed selection
- Stable handling: secure four-wheel chassis for smooth everyday travel
- Made to suit you: personalise with accessories and colour accents

Personalisation

with all interchangeable accent options included as standard, you can mix and match to suit your style and season.



eFoldi Yellow



Moonlight Silver



Harbour Blue



Wherever life takes you, we'll be right here, making sure you have the mobility and freedom to enjoy every moment.

As a recognised leader in the design and manufacture of portable scooters and powered wheelchairs, eFOLDi strives to help people rediscover their freedom and independence.

eFOLDi's award-winning mobility solutions are not only easy to store and transport, they help make travel simple and achievable – whether that's by air, land or sea.

Innovation

Our products are manufactured using the highest quality materials to ensure they're robust enough to withstand the toughest of conditions and terrain. Our frames are made from high-grade aluminium or magnesium alloys, our high-specification brushless motors deliver exceptional performance, and our batteries use advanced lithium-ion technology, so they're not only durable, they offer outstanding levels of reliability.



Cliq Specifications*

Maximum Speed	4mph / 6kph
Maximum Range	9miles / 15km
Maximum User Weight	19 Stone / 120kg / 265lbs
Climbing Ability	1:10 Gradient
Turning Radius	121cm (47.6")
Unfolded Dimensions	100 x 61 x 88cm 39 x 24 x 34.6 inches
Individual Part Weight	3.6kg/7.9lbs 9.3kg/20.5lbs 9.8kg/21.6lbs 7.3kg/15.9lbs
Combined Weight (excl. battery)	28kg/61lbs
Combined Weight (incl. battery)	31kg/69lbs
Type of Use	Pavement, Indoor
Braking System	Electro Magnetic
Airline friendly	Yes

*Specifications may vary due to manufacturing tolerances

7 Year Battery + 2 Year Scooter Warranty***

***terms and conditions apply

Why eFOLDi stands apart




Lightweight
Carry and lift with ease




Warranty
Backed by full support & warranty



Compact
Portable for greater independence



Robust
Peace of mind with strong, durable materials



Scan for more information

eFOLDi

For more information...

Call us today on **0800 999 8123** Or email us at: **sales@efoldi.com**



SunTech UK Limited t/a eFOLDi, Unit H The Courtyard Business Centre, Lonesome Lane, Reigate, Surrey, RH2 7QT

eFOLDi

eFOLDi Cliq

Lightweight mobility solutions

For wherever life takes you

www.efoldi.com



Scan for more information