

Exploded View Compass



Part	SKU	Description	Qty
1	SUN0476	Compass - Throttle	1
2	SUN0274	Lite - Handle bar soft grip	2
3	SUN0484	Compass - Handle bar assembly	1
4	SUN0386	Lite - Tiller Housing with charging port (D-shaped)	1
5	SUN0485	Compass - Forward/Reverse switch	1
6	SUN0486	Compass - Seat support frame welding assembly	1
7	SUN0487	Compass - Steering stem folding lock set with all fittings	1
8	SUN0488	Compass - Stem lower assembly	1
9	SUN0475	Compass- controller	1
10	SUN0491	Compass - front solid wheel 7" cover black	2
11	SUN0490	Compass - Front solid wheel 7"	2
12	SUN0492	Compass -Turning control assembly	1
13	SUN0493	Compass - Backrest pad	1
14	SUN0494	Compass- seat cushion	1
15	SUN0495	Compass - Flexi backsupport	1
16 *	SUN0470	Compass armrest bundle	1
17	SUN0496	Compass - Frame release toggle and nut (yellow)	1
18	SUN0497	Compass - Seat frame support rod	2
19	SUN0069	10ah battery	1
20-1	SUN0499	Compass - Front footrest plastic cover	1
20-2	SUN0498	Compass - Rear footrest plastic cover	1
20-3	SUN0504	Compass - Middle footrest plastic cover	1
21	SUN0500	Compass - Main frame	1
22	SUN0501	Compass - 2" anti-tipping wheel	2
23	SUN0502	Compass - Rear solid wheel 8"	2
24	SUN0503	Compass - Rear solid wheel 8" cover black	2
25	SUN0474	Compass - Motor	1
	* SUN0472	Left armrest	
	* SUN0471	Right armrest	
Not shown:	SUN0477	Compass - main cable harness	