



eFOLDi[®] Information for Airlines and Handlers

Device information

Manufacturer:	eFOLDi
Model:	Cliq
Serial Number:	
Device Value:	£2,299
Maximum Weight:	31kg
Weight Unloaded:	31kg
Height:	88cm
Width:	61cm
Length:	100cm

Manoeuvrability & adjustment features

- ☐ Driving with controller
- ☒ Push in freewheel only
- ☒ Back folds down
- ☐ Back reclines
- ☐ Seat tilts
- ☐ Release lever marked

Battery Information

Battery type:	Lithium-ion
Quantity:	one
Watt-hour (Wh):	288 Wh



eFOLDi[®] Information for Airlines and Handlers

Removeable parts

Component	Remain on device	Cabin
Seat Cushion	✓	
Head Support	N/A	
Arm Supports	✓	
Leg Supports	N/A	
Control Device	✓	
Secondary Control	N/A	
Backrest	N/A	
Tray	N/A	
Belt/Straps	N/A	
Other	N/A	

Passenger Contact

Name: _____

Email: _____

Mobile: _____

Flight Number: _____

CRITICAL: Lift only from base/frame – NEVER from armrests